# TB-IC1020

## Termometro

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### Uso previsto

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Class 1, according to risk.

### Authorization

Designed under Quality Management System standards ISO 13485:2003/NS-EN ISO 13485:2012.

### Characteristics

Dimensions: 70 mm (probe).

Range: -50 to +150°C / -58 to +302°F.

Resolution: 0.1 °C between -19.9 °C to 199.9 °C, otherwise 1 °C.

Sampling Rate: 1 second.

Battery: 1 pc x 1.5 volt type LR44 or equivalent.

Thermometer size: 125 mm (L).

Ambient temperature: -10 °C to 50 °C / +14 °F to +122 °F.

Calibration certificate is provided.

Calibration method: 4 point comparison of temperature reading against reference thermometer.

### Environmental conditions during manufacture

Room temperature 15-35 °C, RH 35-85 %

**Storage conditions**

T = 10-30 °C, RH 30-80 %.

**Transportation conditions**

Storage conditions should be strictly followed. The product should be transported in closed and reinforced boxes in order to avoid damages. Product transportation does not represent any risk for human health.

**Packing**

One unit per box.

Packing information: product code and description, specifications, directions for use, calibration certificate, batch number, serial number and manufacturer information.

**Labelling**

On product: serial number and manufacturer information.

**Possible target markets**

Healthcare, Food, Pharmaceutical and medical devices industries.

**Other important information (optional)**

Use as an accessory for incubators IC10/20, IC10/20FR, IC10/20FRLCD, MiniPro and MiniBio to monitor temperature at  $37 \pm 2$  °C or  $60 \pm 2$  °C.

The metal probe **MUST BE** completely inserted into the incubator.

Note: When necessary, limits of each product's parameters should be specified.

**Description of identified and/or predictable risks which could represent a risk for the product**

Do not bend the probe.

Do not wet.

Do not store at direct sunlight, high temperature or high humidity.

Do not use inside an oven.